

ULTRASOUND MAGIC 3000 PLUS

Mobile & functional



Ultrasound

MOBILE & FUNCTIONAL MAGIC 3000 PLUS

We observe the work of veterinarians in practice and design devices that adapt to the needs of our customers. The MAGIC 3000 PLUS impresses not only with its new design, but also with many new functions and features. Power Doppler (PDI) and PW Doppler allow blood flow to be assessed under a variety of clinical conditions, and the enhanced imaging provides outstanding imaging results for this device class.



Increased picture quality

Modern and user-friendly applications help you to set the optimal picture.

Easy archive

You can easily transfer your pictures, clips, and reports to the practice computer and easily retrieve them from your device via iStation™.

More time for your Patient

The high level of comfort in use and image optimizations save valuable time for your patients.

The MAGIC 3000 PLUS with his compact construction is optimal for use in clinic and my mobile practice.

5,6 kg only



ULTRASOUND MAGIC 3000 PLUS

Areas of Application

- Internal medicine, Gynecology
- Echocardiography

Adjustments for optimal working

- iStation™: intuitive Patient management
- iClear™: real-time image optimization
- iZoom™: automatic zoom
- iTouch™: touch-sensitive controls
- iBeam™: automatic beam steering
- 8-special TGC: 8 special TGC curves
- PW-, Power-Doppler



Technical Data

- LCD-Monitor, ca. 30 cm Diagonal, tiltable
- Image export to Windowsformats
- Intern Storage (320 GB / corresponds to ca. 500.000 pic.)
- Long battery life of up to 2 hours (optional)
- 2 probe-connections
- 2 USB-connections
- Dimensions (in cm):
W 29 x H 35 x D 17

304045

Ultrasound probes

FOR MAGIC 3000 PLUS

Technologie with an impressionnant image quality.

Our transducers are available in several versions, so we can always offer you a solution that meets the needs of your field..



ULTR.PROBE MICRO-CONVEX

Frequency Switching of: 5,0 / 6,5 / 7,5 / 8,5 MHz

303710



ULTR.PROBE LINEAR T-SHAPED

Frequency Switching of: 5,0 / 7,5/8,5/10,0 MHz

303711



ULTR.PROBE LINEAR RECTAL

Frequency Switching of: 4,0 / 5,0 / 6,0 / 7,5 MHz

304021



ULTR.PROBE LINEAR RECTAL

Frequency Switching of: 5,0 / 7,5 / 8,5 / 10,0 MHz

304022

CARRY CASE FOR MAGIC 3000 PLUS

Für den secure transport in ambulant work

304043



TROLLEY FOR MAGIC 3000 PLUS

Easy movement in your treatment room and veterinary clinic.

▶ Printer option

304042

BATTERY FOR MAGIC 3000 PLUS

For 120 Minute continuous work without an electricity source.

304041



Gleichzeitige schwarz-weiß und Power-Doppler Schallung