

ENDOSCOPY EICKVIEW COMPACT SET

Mobile, simple & smart

**Total
control**

*Impressive
versatility*



Endoscopy

MOBILE, SIMPLE & SMART EICKVIEW COMPACT SET

Our devices are designed with the needs of veterinary surgeons in mind. Due to its small diameter and integrated working channel, the specially developed EickView Compact Endoscope provides the perfect solution for rhinoscopy, tracheoscopy and vaginoscopy in small animals.



Intelligent Fusion

Fusion of endoscope and sheath creates a homogenous and minimal invasive distal end.

Clear View

Minimal dirt at distal end.

Easy access

Small diameter enables easy access to most orifices.

High Mobility

Independent operating principle through attachable LED battery light source.

Total control

Integrated working channel and filigree foreign body forceps enable rapid foreign body extractions.

“The EickView Compact Set provides veterinarians with everything that is needed to perform an easy, fast and safe foreign body extraction.”

innovative



ENDOSCOPY

EICKVIEW COMPACT SET

- ▶ Outer Ø: 3.8 mm
- ▶ Working length: 14 cm
- ▶ Working channel Ø: 1.3 mm
- ▶ View angle: 30°
- ▶ Field of view: 90° (wide angle)
- ▶ Endoscope soakable
- ▶ Battery operation time: 4 h
- ▶ Forceps Ø: 1.0 mm
- ▶ Forceps length: 25 cm



- ▶ Scope of delivery: EickView Compact Endoscope incl. otoscope specula kit, LED battery light source, foreign body forceps, charger, spare battery, transport case, two cleaning brushes
- ▶ Stopcocks: Two integrated stopcocks allow steady insufflation and irrigation
- ▶ Light source: Three brightness levels
- ▶ Connection: Storz, Wolf or ACMI connectors
- ▶ Eyepiece: DIN ocular allows the connection of all standard C-mount cameras
- ▶ Padded transport case: W 29 x D 25 x H 9 cm

306122

Endoscopy

ACCESSORIES FOR EICKVIEW COMPACT SET

The following accessories compliment the EickView Compact Set. All listed accessories are of the highest quality and have been tested for longevity.



EICKCAM IV USB HD CAMERA

- For live image reproduction and archiving of images and videos

306182



BIOPSY FORCEPS FOR EICKVIEW

- For extraction of biopsies
- L 25 cm, Ø 1.0 mm

30611804



DISINFECTANT SOLUTION

- Gigasept FF for rigid and flexible endoscopes

307315



DISINFECTANT BASIN

- Disinfectant basin with inner tray, clear lid, 1 l

307365