



EickTron 160 W
Electrosurgical Unit

Electric surgery made easy

EickTron 160 W ElectroSurgical Unit

The EickTron ElectroSurgery Unit is a versatile tool in veterinary medicine allowing it to be used for both monopolar and bipolar instruments. The unit's power adjustment is actuated via a foot or hand switch.

Includes

- ▶ Handpiece with finger switch and connecting cable
- ▶ Electrode set (10 pieces)
- ▶ Cable for neutral electrode
- ▶ Neutral electrode
- ▶ Footswitch
- ▶ Optional Bipolar accessories available

Specifications


- ▶ Input: 220 – 240 V; 50 – 60 Hz
- ▶ Power consumption: 300 VA
- ▶ Frequency: 600 kHz
- ▶ Max. Output 160 W (cut)
- ▶ Max. Output 100 W (coagulate)
- ▶ Max. Output 60 W (bipolar)
- ▶ Includes: Rectal probe
- ▶ Dimensions (in cm): H 10,4 x B 25,4 x T 28,8

323135





EickTron 160 W
Electrosurgical Unit



Accessories for EickTron ElectroSurgical Unit		
Item No.	Description	Pieces
Bipolar Accessories		
323105	EickTron Bipolar Forceps • Autoclavable • Angled • Length: 150 mm	1
323106	EickTron Bipolar Forceps Cable • Autoclavable	1
323139	EickTron Bipolar Adaptor • Autoclavable	1
Monopolar Accessories		
323137	EickTron Monopolar Handpiece • Spare part	1
32313606	EickTron Electrode Angled Hook (L5)	5
32313603	EickTron Electrode (L2)	5
32313607	EickTron Electrode Needle – Angled (L6)	5
32313602	EickTron Electrode Ball – Angled (L10) • Ø 3,0 mm	5
32313608	EickTron Electrode Drop-shaped (L7)	5
32313605	EickTron Electrode Wire Loop (L4) • Ø 8,0 mm	5
32313609	EickTron Electrode Noose (L8)	5
32313604	EickTron Electrode Wire Loop (L3) • Ø 4,0 mm	5
32313601	EickTron Electrode (L1)	5
32313610	EickTron Electrode Ball – Straight (L9) • Ø 3,0 mm	5
323138	EickTron Neutral Electrode Cable • For monopolar electrodes	1
323140	Handle with Finger Switch	1
323141	EickTron Neutral Electrode	1
323143	Monopolar Scissors • Length: 18 cm	1
323144	Connection Cable • For Monopolar Scissors (323143) • Length: 3 m	1