



InfusoVet Easy II
Infusion Pump

veterinary technology for life
Eickemeyer®

In my veins



402080

InfusoVet Easy II Infusion Pump

The EICKEMEYER® InfusoVet Easy II is a programmable infusion pump that combines safety, performance and reliability.

- ▶ Very good value for money
- ▶ LCD display with backlight
- ▶ Intuitive operation with colored, clearly arranged buttons and quick guide on the screen
- ▶ User-friendly through simple programming
- ▶ High patient safety through audible and visual alarms (in case of air bubble, blockage, empty battery, end of infusion)
- ▶ Embedded handle for easy transport
- ▶ The InfusoVet Easy II Infusion Pump (402080) should only be used with EICKEMEYER® Infusion lines (402071) or other pressure/gravity driven infusion line that meets ISO 8536-8 P requirements.

Technical specifications

- ▶ Flow rate: 0,1 up to 1.500 ml/h
- ▶ Accumulated volume: 1 – 9.999 ml
- ▶ Bolus rate from 0,1 to 1.500 ml/h
- ▶ Mains and battery operation, battery life up to 4 hours
- ▶ Alarm volume and alarm sound adjustable
- ▶ With clamp for IV pole
- ▶ Dimensions (in cm): W 14,6 x H 19,3 x D 9,5
- ▶ Weight: 1,6 kg

402080

InfusoVet Cagemount

402076

EICKEMEYER®
Eltastraße 8
78532 Tuttlingen
Germany

T +49 7461 96580 0
F +49 7461 96580 90
E export@eickemeyer.com

TELEPHONE +49 7461 96580 0
www.eickemeyer.com