



EICKEMEYER® Doppler

veterinary technology for life
Eickemeyer®

Accurate & reliable

EICKEMEYER® Doppler

Reliable and accurate non-invasive blood pressure monitor for small animals.

- ▶ Compact, hand-held design
- ▶ High sensibility for a clear signal
- ▶ Adjustable sound volume and headphone connection
- ▶ Focused 8 MHz probe with wide beam for better flow detection in the artery
- ▶ Probe adapted for positioning on the carpal joint or the underside of the tail
- ▶ Rechargeable battery powered, operating time: ± 6 hours
- ▶ Dimensions (in mm): L 175 x W 80 x H 25
- ▶ Lightweight (300 g)

217210



Sphygmomanometer

- ▶ Includes two cuffs: 5.0 and 7.5 cm

217309

Doppler Blood Pressure Cuff, 1.9 cm

- ▶ Ideal for small cats

217307

Doppler Blood Pressure Cuff, 3.3 cm

- ▶ Ideal for cats and small dogs

217308

Doppler Blood Pressure Cuff, 10 cm

- ▶ Ideal for large dogs

217311



217308



217309

EICKEMEYER®
Eltastrasse 8
78532 Tuttlingen
Germany

T +49 7461 96580 0
F +49 7461 96580 91
E export@eickemeyer.com

TELEPHONE +49 7461 96580 0
www.eickemeyer.com