

X-DR Static Classic
Fully digital complete X-ray system



X-DR Static Classic

Complete system

The X-DR Static Classic combines all components and features for high resolution digital images.

The X-ray system and the detector are perfectly matched. The central control system of all components with multi-touch monitor ensures maximum operating comfort.

Powerful X-ray system

Generator: 32kW / max. 450 mA

Voltage range: 40 – 125 kV

Mains connection: 230 V – 50 / 60 Hz, 1-phase


Capacitor technology for higher performance

Excellent X-ray images

Detector 43 cm x 43 cm (17" x 17")

Advanced CsI– scintillator technology ensures brilliant X-ray images.

State-of-the-art image processing optimized for small animals.



Maximum user comfort

Focus fully on the animal

With X-DR Static Classic solution you can focus on your patient.

Optimized hardware

Multi-touch monitor: Zooming or rotating of X-ray images by simply spreading your fingers.

Floating tabletop for simple positioning of the animal before shooting.

Central ON/OFF switch for fast system start.

Intuitive software control

Animal specific organ programs for setting of X-ray system and detector – all settings with just one click.

X-ray images are displayed within few seconds.

Integrated X-AQS viewer facilitates diagnostic viewing directly on the system. Measurement functions included.



Image archive and image transfer

Easy and fast communication with animal owners and colleagues. All necessary functions are integrated:

- Creation of patient CDs
- Image storage on USB stick
- Extensive print function
- Dropbox cloud function optional

Archiving options:

- Integrated image storage – access to archived X-ray images directly on the X-ray system with an option for additional viewer workstations
- Connection to existing image archives



Benefit from many years of experience

- Extensive support with project design
- Modular system configuration
- In-house training for technicians
- Installation support if required
- Second level support



EXAMION GmbH

EXAMION GmbH currently serves an installed base of over 5000 systems. EXAMION is represented by 5 subsidiaries, 46 channel partners and distributors in over 35 countries worldwide. We favor direct contact with our customers, by being available through our extensive network of regional support points and system partners.

Over a network of certified sales and service partners EXAMION provides customers around the world with its latest state-of-the-art X-ray solutions. EXAMION distinguishes itself by its extensive market experience, its advanced cutting-edge technology and the large number of project references in the private and public sector.

EXAMION GmbH (Headquarter)

Erich-Herion-Str. 37 Telefon +49 711 1200020 www.examion.com
70736 Fellbach, Germany Fax +49 711 12000222 info@examion.com

Certified Partner of EXAMION



Société NOVEDA
Office: 8 Rue Bouhajla, Ariana 2080 - Tunisia
Sales: Avenue de l'environnement, Cité Essaada Nr. 157,
Sidi Thabet 2020 - Tunisia