




Veterinary Pulse  
Oximeter LifeVet P

veterinary technology for life  
**Eickemeyer**®

# Every breath you take.



321865



321866

## Veterinary Pulse Oximeter LifeVet P

The LifeVet P measures and displays oxygen saturation (SpO<sub>2</sub>) and pulse rate. An LCD screen displays the numeric values and waveforms from these parameters. Suitable for intermittent or continuous monitoring of canine, feline and equine cases.

- ▶ Displays SpO<sub>2</sub> number and waveform simultaneously
- ▶ Long battery life – up to 48 hours on 4xAA NiMH rechargeable batteries
- ▶ Automatic power-off function for power saving
- ▶ Compact, lightweight & easy to clean
- ▶ No calibration required
- ▶ Data can be transferred to PC for storage, review and printing (Optional part required)
- ▶ Large data storage capacity
- ▶ Audible and visible alarms
- ▶ Displays saturation level: 1 % – 100 %
- ▶ SpO<sub>2</sub> resolution: 1 %
- ▶ Displays pulse rate: 25 – 350 bpm
- ▶ Pulse rate resolution: 1 bpm
- ▶ Includes two sensor clips (small and large)
- ▶ Includes charging stand and Ni-MH rechargeable AA battery pack

### 321865

Pulse Oximeter LifeVet P Protective Cover

- ▶ Colour: green

### 321866

**EICKEMEYER**®  
Eltastrasse 8  
78532 Tuttlingen  
Germany

T +49 7461 96580 0  
F +49 7461 96580 91  
E [export@eickemeyer.com](mailto:export@eickemeyer.com)

**TELEPHONE +49 7461 96580 0**  
[www.eickemeyer.com](http://www.eickemeyer.com)