


ANAESTHESIA LifeVet Monitors


Accurate and portable



Anaesthesia

ACCURATE AND PORTABLE LifeVet M / LifeVet C

Our devices are designed with the needs of veterinary surgeons in mind. LifeVet monitors are helpful for straightforward and reliable multi-parameter assessment during anaesthesia.



Rapid Documentation

The included software allows paperless storage of data.

Information Display

The high resolution screen allows clear visualisation of patients physiological parameters to provide continuous real time data.

Enhanced patient safety

The optional ETCO₂ measurement feature via the side stream allows even greater levels of monitoring during anaesthesia.

Compact Design

The machine's innovative design allows you to monitor patients in restricted spaces.

"With the LifeVet monitors I have peace of mind when my patients are under anaesthesia."


Dependable



ANAESTHESIA

LifeVet M/LifeVet C

- ▶ Pulse Oximetry with adjustable alarm
- ▶ Accurate pulse rate and SpO₂ display
- ▶ Bar graph to display pulse intensity
- ▶ 5 Lead ECG for monitoring cardiac function
- ▶ Oral or rectal temperature measurement
- ▶ Noninvasive blood pressure monitoring
- ▶ Battery or mains operated
- ▶ Later upgrade to capnography possible



321872



321870

- ▶ 10.1" High resolution colour display
- ▶ SpO₂ Sensors (Small and Large)
- ▶ Documentation software included
- ▶ Includes CO₂ preparation for capnography upgrade
- ▶ Includes alligator clips
- ▶ Includes ECG cable, NiBP cable and extension cable
- ▶ Includes blood pressure cuffs (cm): 3-6/6-11/8-15 (W 2.4/4.0/5.7)
- ▶ Includes Lithium ion battery
- ▶ Dimensions: 31.8 x 26.5 x 13.7 cm (W x H x D)
- ▶ Weight: approx. 4.2kg

321870

- ▶ LifeVet M includes CO₂ module, CO₂ moisture measurement tube and a connector adaptor

321872

Anaesthesia

ACCESSORIES FOR LifeVet Monitors

The following accessories support the use of LifeVet monitors. All listed accessories are of the highest quality and have been tested for longevity.



PULSE OXIMETRY SpO₂ SENSOR

- Includes a small and large sensor to accommodate different animal sizes


321868



LifeVet C (CO₂ MODULE UPGRADE)

- ETCO₂ sidestream mode helps monitor respiratory rate
- Plethysmograph sensor for ETCO₂ (capnography)

321878



LifeVet MONITOR MOBILE STAND

- Stand and tray constructed from aluminium
- Basket constructed from powder coated iron

321873



WALL MOUNT FOR LifeVet MONITOR

- With Basket

321875



LifeVet MONITOR BLOOD PRESSURE CUFFS

- Circumference 3 – 6 cm
E32187003
- Circumference 4 – 8 cm
E32187013
- Circumference 6 – 11 cm
E32187002
- Circumference 7 – 13 cm
E32187014
- Circumference 8 – 15 cm
E32187004

CAPNOGRAPH CONNECTOR FOR CO₂-MODULE

- Connects endotracheal tube to Y-connector
- 2 per pack
321373

FILTER FOR CO₂-MODULE

- Damp- and particle filter
E32187803

DRYING LINE FOR CO₂-MODULE

- Always use with CO₂ Module
E32187802

CO₂ SAMPLE LINE FOR CO₂-MODULE

- With Luer-Lock connections
E32187801