

## Oxygen Therapy Kit



A 'must have' for emergency service provider and lone worker, a veterinary professional may provide oxygen therapy single-handed.

### Oxygen Therapy Box

A simple and effective way to provide oxygen or nebulising therapy, this is a cleverly designed PVC canopy which encloses a basket and attaches to an oxygen source.

- Instant oxygen or nebulising therapy
- Full patient visibility through the 3 mm PVC canopy for observation of respiratory patterns during treatment
- Economical use of oxygen

The welded PVC seams keep the interior airtight for initial oxygenation. Venting access via the zippered top (and handle opening) allows the release of excess oxygen and exhaled carbon dioxide. Venting can be altered depending on the patients' condition and oxygen percentage requirements. Remembering that most respiratory patients do much better without restraint (stress plays a major roll in exacerbating symptoms) the animal may be monitored without the anxiety associated with a nasal cannula or masking.

The package includes one used, reconditioned OxyVet II Oxygen Concentrator.

### OxyVet Oxygen Concentrator

Transfers ambient air into medical oxygen, quiet operation as a whisper at only 40 db(A)

- Up to 95 % pure oxygen, 5 l/min, incl. flowmeter and humidifier for oxygen therapy
- Highest patient security by visual and acoustic alarm functions
- Comfortable, time-saving and cost-effective
- Designed for continuous operation. Not suitable for use with anesthesia devices.

#### EICKEMEYER®

Eltastraße 8  
78532 Tuttlingen  
Germany

T +49 7461 96580 0  
F +49 7461 96580 90  
E [export@eickemeyer.com](mailto:export@eickemeyer.com)

TELEPHONE +49 7461 96580 0

[www.eickemeyer.com](http://www.eickemeyer.com)